FertiCult™ Flushing medium EN

Cell culture medium for washing of human ova, spermatozoa and embryos, for swim-up of spermatozoa, sperm injection in ICSI, intra uterine insemination and embryo transfer

<p>Doc. reference: FP09 I08 R01 C.4 Update: 24.10.2014</p>
<p>USED ABBREVIATIONS</p>
<p>ICSI Intracytoplasmatic Sperm Injection IUI Intra uterine insemination</p>
<p>GENERAL INFORMATION AND INTENDED USE</p>
<p>FertiCult Flushing medium is a formulation for washing of human ova, spermatozoa and embryos. FertiCult Flushing medium can also be used for swim-up techniques of human spermatozoa, sperm injection in oocytes during ICSI, introduction of washed spermatozoa in the uterus (IUI) and for embryo transfer. The medium is complete and needs no further additives. The medium contains HEPES; no CO₂ incubation is required. FertiCult Flushing medium does not contain heparin, thus for oocyte pick-up, use FertiCult Aspiration medium. As with all IVF media, FertiCult Flushing medium ought to be preincubated in the incubator for 12 hours before use (with lid closed).</p>
<p>COMPOSITION</p>
<p>FertiCult Flushing medium is a ready-to-use HEPES-buffered medium which also contains bicarbonate, physiologic salts, glucose, lactate and human serum albumin (4.00g/liter). Gentamicin Sulphate can be added upon request (10mg/liter).</p>
<p>MATERIAL INCLUDED WITH THE KIT</p>
<p>FLUSH020 » 5 x 20ml of FertiCult Flushing medium FLUSH22.8 » 1 x 22.8ml of FertiCult Flushing medium FLUSH050 » 5 x 50ml of FertiCult Flushing medium FLUSH100 » 3 x 100ml of FertiCult Flushing medium FLUSH_002_PHR » 1 x 2ml of FertiCult Flushing medium with Phenol Red FLUSH020PHR » 1 x 2ml of FertiCult Flushing medium with Phenol Red FLUSH050PHR » 5 x 50ml of FertiCult Flushing medium with Phenol Red FLUSH100PHR » 3 x 100ml of FertiCult Flushing medium with Phenol Red FLUSH500PHR » 1 x 500ml of FertiCult Flushing medium with Phenol Red FLUSH020PRG » 5 x 20ml of FertiCult Flushing medium with Phenol Red and Gentamicin FLUSH050PRG » 5 x 50ml of FertiCult Flushing medium with Phenol Red and Gentamicin FLUSH100PHR_G » 3 x 100ml of FertiCult Flushing medium with Phenol Red and Gentamicin FLUSH500PHR_G » 1 x 500ml of FertiCult Flushing medium with Phenol Red and Gentamicin WASH005 » 1 x 5ml of Washing/Insemination medium WASH020 » 1 x 20ml of Washing/Insemination medium</p>
<p>MATERIAL NOT INCLUDED WITH THE KIT</p>
<p>» Incubator at 37°C (No CO₂) » Petri dishes » Microscope » Test tubes » LAF bench (ISO Class 5) » Syringe (e.g. 1ml Plastipack) » Catheter</p>
<p>PRODUCT SPECIFICATIONS</p>
<p>» Chemical composition » pH: between 7.30-7.60 (Release criteria) » Osmolality: 270-290 mOsm/kg » Sterility: sterile (SAL 10⁻³) » Endotoxins: < 0.25 EU/ml » Mouse Embryo Assay (blastocysts after 96h): ≥ 80% after 30min of exposure (zygote stage) » Use of Ph Eur or USP grade products if applicable » Certificate of analysis and MSDS are available upon request</p>

<p>PRE-USE CHECKS</p>
<p>» Do not use the product if it becomes discoloured (if medium contains phenol red), cloudy, or shows any evidence of microbial contamination. » Do not use the product if seal of the container is opened or defect when the product is delivered.</p>
<p>STORAGE INSTRUCTIONS</p>
<p>» Store FertiCult Flushing medium without gentamicin between 2-25°C. » Store FertiCult Flushing medium with gentamicin between 2-8°C. » Once opened store all products between 2-8°C. » Do not freeze before use. » Keep away from (sun)light. » The products can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C. » Do not use after expiry date. » Stable after transport (max. 5 days) at elevated temperature (≤ 37°C).</p>
<p>WARNINGS AND PRECAUTIONS</p>
<p>Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens. Always work under strict hygienic conditions (e.g. LAF-bench ISO Class 5) to avoid possible contamination, even when FertiCult Flushing medium contains Gentamicin. FertiCult Flushing with phenol red and gentamicin contains the antibiotic Gentamicin Sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.</p>
<p>METHOD</p>
<p>Washing of spermatozoa (suggested procedure)</p>
<p>The washing of spermatozoa can be done at room temperature or at 37°C.</p> <ol style="list-style-type: none">Add 5ml FertiCult Flushing medium to the native semen sample and mix. Centrifuge for 15 minutes at approximately 300g. Remove supernatant and leave about 0.5ml of semen in the centrifuge tube. Add 5ml FertiCult Flushing medium to the test-tube. Mix the solution gently until the pellet is completely dissolved. Centrifuge again for 10 minutes at 300g. Perform a swim-up.
<p>Swim-up procedure (according to WHO, 2010)</p>
<ol style="list-style-type: none">Gently layer 1,5ml FertiCult Flushing medium over 1ml of washed semen in a conical based centrifuge tube. Incline the tube at an angle of 45° and incubate for 1 hour at 37°C. Gently turn to the upright position and remove the uppermost 1ml. Dilute this aliquot of motile cells with 8 volumes FertiCult Flushing medium. Centrifuge for 15 minutes at 300g, remove supernatant and finally resuspend in 0.5ml of FertiCult Flushing medium.
<p>Embryo transfer (according to Brinsden, 2005)</p>
<ol style="list-style-type: none">Take a sterile syringe, fill with FertiCult Flushing medium, and eject air bubbles. Attach the syringe to a catheter, and eject the medium. Draw up FertiCult Flushing medium into the syringe and push down the piston to the 10µl calibration mark. Draw up the embryo(s) into the catheter so that the volume to be transferred is approximately 20-30µl. Hand catheter and syringe to the clinician for insertion and transfer in the uterus.

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- WHO laboratory manual for the examination of human semen and sperm-cervical mucus interaction, 3rd edition, Cambridge University Press, Avon, 1992, pp. 107.
- WHO Laboratory Manual for the Examination of Human Semen and semen-cervical mucus interactions, Cambridge University Press, 1999.
- Brinsden P.R. (2005). Bourn Hall Clinic Procedures and Protocols. A textbook of in vitro fertilization and assisted reproduction, 3rd edition, Taylor & Francis. Procedure No. K.7.

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