SpermMar Test IgG

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1. BOETTCHER B., HJORT T., RUMKE Ph., SHULMAN S. and VYAZOV O.E. (eds.): Auto and iso-antibodies to antigens of the human reproductive system. 1. Results of an international comparative study of antibodies to spermatozoa and other antigens detected in sera from infertile patients deposited in the WHO Reference bank for reproductive immunology. Acta Pathol Microbiol Scand, 1977, 258 (suppl.) : 1-69
2. COMHAIRE F.H. and KUNNEN M.: Factors affecting the probability of conception after treatment of subfertile men with varicocele by transcatheter embolisation with Bucrylate. Fertil Steril, 1985, 43: 781-786.
3. FRIBERG J.: Immunoglobulin concentration in serum and seminal fluid from men with and without sperm-agglutinating antibodies. Obstet. Gynecol., 1980, 136 : 671-675.
4. GOLOMB J., VARDINON H., HOMONNAI Z.T., BRAF Z. and YUST I.: Demonstration of antispermatozoal antibodies in varicocele-related infertility with an enzymelinked Immunosorbent assay (ELISA). Fertil Steril, 1986, 45 : 397-402.
5. JAGER S., KREMER J. and VAN SLOCHTEREN-DRAAISMA T.: A simple method of screening for antisperm antibodies in the human male. Detection of spermatozoal surface IgG with the direct mixed antiglobulin reaction carried out on untreated fresh human semen. Int. J. Fertil., 1978, 23 : 12-21.
6. JAGER S., KREMER J., KUIKEN J. and VAN SLOCHTEREN-DRAAISMA T.: Immunoglobulin class of antispermatozoal antibodies from infertile men and inhibition of in vitro sperm penetration into cervical mucus, IntJ. Androl., 1980, 3: 1-14.
7. RUMKE P.: The origin of immunoglobulins in semen. Clin. Exp. Immunol., 1976, 17 : 287-297.
8. SHULMAN J.F. and SHULMAN S.: Methylprednisolone treatment of immunologic infertility in the male. Fertil Steril, 1982, 38 : 591-599.
9. STEDRONSKA J. and HENDRY W.F.: The value of the mixed antiglobulin reaction (MAR-Test) as an addition to routine seminal analysis in the evaluation of the subfertile couple. Am. J. Reprod. Immunol., 1983, 3 : 89-91.
10. VERMEULEN L. and COMHAIRE F.H.: Le test <MAR> aux particules de Latex, et le test spermatoxique selon Suominen : simplification et nouveautй dans l'arsenal du diagnostique immunologique. Contraception-Fertilite-Sexualite, 1983, 11 (suppl.) : 381-384.

11. W.H.O. 1984: Workshop on the standardized investigation of the infertile couple, moderator P. Rowe, coordinator M. Darling. Fertil Steril (eds) R.F. HARRISON, J. BONNAR and W. THOMPSON. Publ. M.T.P.-Press Ltd. (Lancaster, Boston, The Hague, Dordrecht), 1984, 427-431.
12. W.H.O.: Laboratory manual for the examination and processing of human semen. W.H.O., 5th Edition, 2010.
13. WORLD HEALTH ORGANIZATION: Towards more objectivity in diagnosis and management of male infertility. Int J Androl, 1987, Suppl 7.
14. KREMER J. and JAGER S.: The significance of antisperm antibodies for sperm cervical mucus interaction. Hum. Reprod., 1992, 7 : 781-784.
15. KAY D.J. and BOETTCHER B.: Comparison of the SpermMar test with currently accepted procedures for detecting human sperm antibodies. Reprod. Fer. Dev., 1992, 4, 175-181.
16. ACKERMAN S., Mc GUIRE G., FULGHAM D.L. and ALEXANDER N.: An evaluation of a commercially available assay for the detection of antisperm antibodies. Fertil Steril 1988, 49 : 732-734.
17. BRONSON R, COOPER G and ROSENFELD D: Sperm antibodies: their role in infertility. Fertil Steril, 1984, 42: 171-183.
18. CLARKE GN, STOJANOFF A, CAUCHI MN and JOHNSTON WIH: The immunoglobulin class of antispermatozoal antibodies in serum. Am J Reprod Immunol Microbiol, 1985, 7: 143-147.
19. HINTING A, VERMEULEN L and COMHAIRE F: The indirect mixed antiglobulin reaction test using a commercially available kit for the detection of antisperm antibodies in serum. Fertil Steril, 1988, 49 :1039-1044.
20. MEINERTZ H and HJORT T: Detection of autoimmunity to sperm: mixed antiglobulin reaction (MAR) test or sperm agglutination? A study on 537 men from infertile couples. Fertil Steril, 1986, 46 : 86-91.
21. PARSLAW JM, POULTON TA, BESSER GM and HENDRY WF: The clinical relevance of classes of immunoglobulins on spermatozoa from infertile and vasectomized males. Fertil Steril, 1985, 43 : 621.
22. RUMKE P and HELLINGA: Autoantibodies against spermatozoa in sterile men. Am J Clin Pathol, 1959, 32 : 357-363.
23. UELING DT: Secretory IgA in seminal fluid. Fertil Steril, 1971, 22 : 769-773.
24. WILSON L: Sperm agglutinins in human semen and blood. Proc Soc Exp Biol Med, 1954, 85 : 652-655.

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1. BOETTCHER B., HJORT T., RUMKE Ph., SHULMAN S. and VYAZOV O.E. (eds.): Auto and iso-antibodies to antigens of the human reproductive system. 1. Results of an international comparative study of antibodies to spermatozoa and other antigens detected in sera from infertile patients deposited in the WHO Reference bank for reproductive immunology. Acta Pathol Microbiol Scand, 1977, 258 (suppl.) : 1-69
2. COMHAIRE F.H. and KUNNEN M.: Factors affecting the probability of conception after treatment of subfertile men with varicocele by transcatheter embolisation with Bucrylate. Fertil Steril, 1985, 43: 781-786.
3. FRIBERG J.: Immunoglobulin concentration in serum and seminal fluid from men with and without sperm-agglutinating antibodies. Obstet. Gynecol., 1980, 136 : 671-675.
4. GOLOMB J., VARDINON H., HOMONNAI Z.T., BRAF Z. and YUST I.: Demonstration of antispermatozoal antibodies in varicocele-related infertility with an enzymelinked Immunosorbent assay (ELISA). Fertil Steril, 1986, 45 : 397-402.

5. JAGER S., KREMER J. and VAN SLOCHTEREN-DRAAISMA T.: A simple method of screening for antisperm antibodies in the human male. Detection of spermatozoal surface IgG with the direct mixed antiglobulin reaction carried out on untreated fresh human semen. *Int. J. Fertil.*, 1978, 23 : 12-21.
6. JAGER S., KREMER J., KUIKEN J. and VAN SLOCHTEREN-DRAAISMA T.: Immunoglobulin class of antispermatozoal antibodies from infertile men and inhibition of in vitro sperm penetration into cervical mucus, *IntJ. Androl.*, 1980, 3: 1-14.
7. RUMKE P.: The origin of immunoglobulins in semen. *Clin. Exp. Immunol.*, 1976, 17 : 287-297.
8. SHULMAN J.F. and SHULMAN S.: Methylprednisolone treatment of immunologic infertility in the male. *Fertil Steril*, 1982, 38 : 591-599.
9. STEDRONSKA J. and HENDRY W.F.: The value of the mixed antiglobulin reaction (MAR-Test) as an addition to routine seminal analysis in the evaluation of the subfertile couple. *Am. J. Reprod. Immunol.*, 1983, 3 : 89-91.
10. VERMEULEN L. and COMHAIRE F.H.: Le test <MAR> aux particules de Latex, et le test spermatoxique selon Suominen : simplification et nouveauté dans l'arsenal du diagnostique immunologique. *Contraception-Fertilite-Sexualite*, 1983, 11 (suppl.) : 381-384.
11. W.H.O. 1984: Workshop on the standardized investigation of the infertile couple, moderator P. Rowe, coordinator M. Darling. *Fertil Steril* (eds) R.F. HARRISON, J. BONNAR and W. THOMPSON. Publ. M.T.P.-Press Ltd. (Lancaster, Boston, The Hague, Dordrecht), 1984, 427-431.
12. W.H.O.: Laboratory manual for the examination and processing of human semen. W.H.O., 5th Edition, 2010.
13. WORLD HEALTH ORGANIZATION: Towards more objectivity in diagnosis and management of male infertility. *Int J Androl*, 1987, Suppl 7.
14. KREMER J. and JAGER S.: The significance of antisperm antibodies for sperm cervical mucus interaction. *Hum. Reprod.*, 1992, 7 : 781-784.
15. KAY D.J. and BOETTCHER B.: Comparison of the SpermMar test with currently accepted procedures for detecting human sperm antibodies. *Reprod. Fer. Dev.*, 1992, 4, 175-181.
16. ACKERMAN S., Mc GUIRE G., FULGHAM D.L. and ALEXANDER N.: An evaluation of a commercially available assay for the detection of antisperm antibodies. *Fertil Steril* 1988, 49 : 732-734.
17. BRONSON R, COOPER G and ROSENFIELD D: Sperm antibodies: their role in infertility. *Fertil Steril*, 1984, 42: 171-183.
18. CLARKE GN, STOJANOFF A, CAUCHI MN and JOHNSTON WIH: The immunoglobulin class of antispermatozoal antibodies in serum. *Am J Reprod Immunol Microbiol*, 1985, 7: 143-147.
19. HINTING A, VERMEULEN L and COMHAIRE F: The indirect mixed antiglobulin reaction test using a commercially available kit for the detection of antisperm antibodies in serum. *Fertil Steril*, 1988, 49 :1039-1044.
20. MEINERTZ H and HJORT T: Detection of autoimmunity to sperm: mixed antiglobulin reaction (MAR) test or sperm agglutination? A study on 537 men from infertile couples. *Fertil Steril*, 1986, 46 : 86-91.
21. PARSLAW JM, POULTON TA, BESSER GM and HENDRY WF: The clinical relevance of classes of immunoglobulins on spermatozoa from infertile and vasectomized males. *Fertil Steril*, 1985, 43 : 621.
22. RUMKE P and HELLINGA: Autoantibodies against spermatozoa in sterile men. *Am J Clin Pathol*, 1959, 32 : 357-363.
23. UELING DT: Secretory IgA in seminal fluid. *Fertil Steril*, 1971, 22 : 769-773.
24. WILSON L: Sperm agglutinins in human semen and blood. *Proc Soc Exp Biol Med*, 1954, 85 : 652-655.

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